


Abstract

ChameLogK is a cutting-edge method for measuring chameleonicity, a vital factor in assessing cell permeability propensity of small molecules, including 'Beyond Rule-of-Five' such as PROTACs, macrocycles and peptides. Chameleonicity is ability of molecules to conformationally adapt to their environment, driven by intramolecular hydrogen bonds; ChameLogK combines lipophilicity (BRlogD) and polarity ($\Delta \log k_w^{IAM}$) descriptors. Leveraging literature methods, o2h has developed a reliable, high-throughput platform to rapidly and inexpensively measure chameleonicity, accelerating optimization of Beyond Rule-of-Five therapeutic compounds

Chameleonicity: An indicator of bioavailability and cell permeability for Beyond Rule of Five molecules

- Chameleonicity is the capacity of a molecule to hide polarity in non-polar environments and expose it in water, help achieving sufficient permeability and solubility for drug molecules with high molecular weight [1]




Adapted from Danielus et al, Nature 2020 [2]

- For many years, it has been known that bRo5 molecules, such as macrocycles and small peptides, can exhibit chameleonic properties – conformational flexibility - driven by the formation of intramolecular hydrogen bonds [3]

- bRo5 compounds often suffer from a poor permeability and bioavailability profile. The design and optimization of compounds with chameleonic properties is central to optimizing solubility and permeability [4]

- BRlogD, $\Delta \log k_w^{IAM}$ and ChameLogK are three complementary measures revealing lipophilicity, solubility, exposed polarity, and chameleonic masking; together forming a triad of drug-likeness, guiding smarter compound prioritization by explaining permeability and bioavailability [5]



- Other established phys chem measures relating to permeability, such as ePSA [6][1] are comparatively low throughput and costly, revealing a need for an easily-accessible alternative measure to support rapid design and optimization of bRo5 compounds

- The ChameLogK method was pioneered in by Giulia Caron's group in Turin [5] and has garnered significant interest in the academic setting, but had not been developed into a fully-transferable platform for routine screening in the commercial setting

- o2h discovery is the first in the world to offer the ChameLogK assay to clients to support their drug discovery programs

Assay Validation and Platform Development

o2h conducted a rigorous validation process, with a variety of different chemical species, underscoring the reliability of the assay as a robust and transferable platform, broadly applicable for multiple drug discovery programs.

- Validated with >20 different compounds, including macrocycles, non-macrocylic compounds, and PROTACs
- >100 compounds examined to date confirm versatility and effectiveness in designing bioavailable compounds beyond Ro5 constraints




Figure 1: A. Output data showing correlation of Log K' with % acetonitrile for a PROTAC and a classic Ro5 compound; B. Summary of chameleonicity across a range of chemical species

Results

- PROTACs:** Displayed distinct "Chameleonic" qualities.
- Non-macrocylics:** Moderate chameleonicity, showing expected flexibility compared to classic Ro5 and PROTACs
- Macrocylics:** Led by Cyclosporin, showed variable behaviour
- E3 Ligands:** Consistent behaviour, possibly involving internal folding or interactions
- Ro5 molecules** like caffeine, and diazepam defied expectation

o2h has streamlined processes for **maximal throughput and productivity**, creating a **cost-effective and easy-access** version of this ground-breaking technique


Measuring ChameLogK: How it works

ChameLogK is the difference ($\Delta \log k'$) between the experimental retention at 100% MeCN on a PLRP-S column, and the extrapolated retention at 100% MeCN obtained from the linear region measured at 50/60/70% MeCN. Positive values indicate extra retention in the non-polar limit, consistent with polarity masking/folding (i.e. chameleonicity).

Method features:

- HPLC-based method, using PLRP-S column
- Record the retention time (RT) of compound in differing ratios of water:acetonitrile
- Linear regression to fit of Log k' PLRP-S and the % acetonitrile; Chameleonic properties are observed as a steeper downward slope with an up-tick
- Extrapolate to calculate at 100% acetonitrile and calculate ChameLogK using the formula:

$$\text{ChameLogK} = \text{Exp. } \log k'_{100} - \text{Ext. } \log k'_{100}$$



Graphics adapted from [2] and [5]


Supporting a range of modalities

- The **increasing need to explore the bRo5 space** to access novel chemistry, better mechanisms of action, and exciting new modalities, requires a way to quickly and inexpensively measure chameleonicity
- This assay represents a **significant step forward** in assessing molecules such as PROTACs and optimization of their properties [7] Enabling:



A key part of your drug discovery program

- o2h offers the ChameLogK assay as part of a comprehensive suite of ADMET and assay biology, enabling full characterisation of compounds



- o2h also offers a range of integrated services linking directly to this field:



- PROTACs:** off-the-shelf-toolbox, featuring >100 building blocks, a range of E3 ligands and ready-to-use linkers
- Custom Peptide Synthesis:** parallel and solid & solution phase; 'peptide guarantee'
- Macrocylics** and natural product chemistry
- Molecular glues**, with comprehensive assay suite for in-cell target engagement and degradation
- Medicinal and synthetic** chemistry

References

- David et al (2021) 'Prediction of Chameleonic Efficiency', ChemMedChem 16(7), pp2669-2685
- Danielus et al, (2020) 'Solution Conformations Explain the Chameleonic Behaviour of Macrocylic Drugs', Nature 26(23), pp5231-5244
- Bockus and Lokey, (2017) 'Chapter 5: Bioactive and Membrane-Permeable Cyclic Peptide Natural Products', Practical Medicinal Chemistry with Macrocylics: Design, Synthesis, and Case Studies, John Wiley and Sons, Inc.
- Ramelot et al, (2023), 'Cell-permeable chameleonic peptides: Exploiting conformational dynamics in de novo cyclic peptide design, Current Opinion in Structural Biology 80, p102603
- Garcia Jimenez et al (2023), 'ChameLogK: A Chromatographic Chameleonicity Quantifier to Design Orally Bioavailable Beyond-Rule-of-5 Drugs', Journal of Medicinal Chemistry 66(15), p10681. **Seminal article by Giulia Caron's group**
- Goetz et al (2014), 'ePSA: A Novel Supercritical Fluid Chromatography Technique Enabling the Design of Permeable Cyclic Peptides', ACS Medicinal Chemistry Letters, 5(10).
- Apprato et al (2024), 'Exploring the chemical space of orally bioavailable PROTACs', Drug Discovery Today 29(4), p103917.

Want to learn more about this assay and how o2h supports integrated drug discovery programs...?

Scan the QR code to explore our website:



Read our blog on ChameLogK:



Contact us: info@o2h.com

